



Material Safety Data Sheet

SECTION I - Product/ Company Identification

Product Trade Name: TRULOCK PRIMER ACTIVATED ADHESIVE
Supplier Name: RMO, INC.
Address: 650 West Colfax Avenue
City State Zip: Denver, CO 80204 U.S.A.
Emergency Telephone Number: (800) 424-9300 (Chemtrec)
Information Telephone Number: (800) 835-0885 / (847) 677-5500
Date Revised: 04-Jan-2016 Rev 001

SECTION II - Hazardous Ingredients/Identity Information

Hazardous Components:	Concentration	CAS NO.	OSHA PEL	ACGIH TLV	OTHER LIMITS
Bis GMA	15 - 25 %	1565-94-2	N/E	N/E	
Silica	65 - 80 %	7631-86-9	N/E	N/E	
Alkyl dimethacrylate	0 - 10 %	109-16-0	N/E	N/E	
Peroxide Initiator	< 5 %	94-36-0	N/E	N/E	

SECTION III - Hazards Identification

HAZARD CLASSIFICATION Harmful if swallowed

SECTION IV - First Aid Measures

Skin: Wash off with soap and water. If irritation persists consult a physician.
Ingestion: Do not induce vomiting, seek medical attention.
Eyes: Flush with plenty of water. Seek medical attention by an Ophthalmologist.
Inhalation: N/A

SECTION V - Fire-fighting Measures.

Suitable extinguishing Media: Chemical Foam, Carbon Dioxide, or Dry Chemical
Forbidden extinguish Media: Unknown
Special Exposure Hazards: None. Wear self contained breathing equipment
Unusual fire and explosion hazards: Heat can initiate exothermic polymerization.
Special Protection Equipment: Unknown

SECTION VI - Accidental release measures

Personal Precautions: See sections VII and VIII
Environmental Precautions: Do not dump in sewer any spilled material.
Methods for Cleaning up: Clean up using impervious gloves and dispose of in an approved manner.



Material Safety Data Sheet

SECTION VII - Handling and Storage

Handling: Avoid unnecessary exposure.
Storage: Room Temperature or lower. Do not freeze, Keep dry and away from sources of heat.
Containers: Supplied by manufacturer.
Incompatibilities: Amines and reducing agents
Other Precautions: Use according to directions. Avoid contact with eyes, skin and soft tissue.

SECTION VIII - Exposure Controls / Personal Protection.

General Measures: Ventilation: good general ventilation
Limits: TWA: not determined; PEEL not determined
Respiratory protection: Generally not necessary
Hands Protection: Impervious Gloves recommended
Eye Protection: Safety Glasses or Goggles recommended
Skin Protection: Respect good personal hygiene practice. Avoid unnecessary exposure. Other
Equipment: Lab Coat.

SECTION IX - Physical and Chemical Properties.

Appearance:	Viscous tinted paste	Explosion Limits:	Not Determined
		Flash Point	Not determined
Odor	Slight acrylic	Water solubility:	Insoluble
Vapor Pressure	Not Determined	Specific Gravity	Not Determined
Boiling Point:	Not Determined	Evaporation Rate:	Not Determined

SECTION X - Stability and Reactivity.

Stability: Stable in normal storage conditions Conditions to
avoid: Temperatures over 40 °C and contamination
Materials to avoid: Amines and reducing agents
Hazardous Decomposition Products: Carbon Oxides
Hazardous Polymerization Products: May occur

SECTION XI - Toxicological Information.

Ingestion: May cause indigestion
Eye contact: May cause irritation, pain, eye burns or permanent eye damage.
Skin Contact: May cause irritation or rash for persons with preexistent skin sensitivity.
Inhalation: Inhalation may cause irritation of the respiratory system.
Prolonged exposure: N/A

SECTION XII - Ecological Information.

There is no specific information on the product. Use according to good working practice. Avoid releasing the product into environment.

SECTION XIII - Disposal Considerations.

Incinerate or dispose according to local regulations.



Material Safety Data Sheet

SECTION XIV - Transport Information.

Substances determining Hazard: Methacrylic Monomers
SEA TRANSPORTATION I
MCO # N/A
UN # N/A
AIR TRANSPORTATION
ICAO/IATA CLASS: N/A
GROUND TRANSPORTATION:
RID/ADR Class N/A
Substance Identification #: N/A
Item # N/A
Hazard Label: N/A
Hazard Identification #: N/A

SECTION XV - Regulatory Information

CLASSIFICATION IN ACCORDANCE WITH EEC DIRECTIVES.
Hazard Classification Code: Xi, Xn Irritant and harmful
Risk Phrases R22 Harmful if swallowed
R36/37/38 Irritating to eyes, respiratory system and skin
Safety Phrases S15 Keep away from heat
S23 Do not breathe vapor
S24/25 Avoid contact with skin and eyes
S26/28 Wash with plenty of water
Ingredients determining the label: Methacrylate Esters Monomers

SECTION XVI - Other Information

None.

CAUTION: PRODUCT FOR PROFESSIONAL USE ONLY.

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.