

FLI® Signature Clear Brackets

FLI® Signature Clear Brackets offer natural aesthetic appeal to patients with the same reliable functionality preferred by doctors.

Ligation: Tie wings incorporate undercuts and downdraft for easy and secure ligation


Finish: Fully transparent for unmatched aesthetics

Reduced Friction: Smooth rounded flared archslot lead-ins for reduced friction

Finish: Smooth rounded surface for exceptional patient comfort

Base: Patented mushroom style base rails deliver superior mechanical bonding

Base: Contoured design provides optimum adaption to tooth surfaces



MATERIAL: Polycrystalline Clear Ceramic – 99.9% Pure Alumina Oxide

PRESCRIPTIONS: Roth, MBT, & Bioprogressive

SLOT SIZE: .018 & .022

FEATURES

BENEFITS

Tie wings with deep undercuts and downdraft	Easy, more secure ligation
Fully transparent	Unmatched aesthetics and true natural appearance
Smooth, polished surfaces	No cloudy appearance as compared to traditional ceramics
Low-profile with rounded contours	Exceptional patient comfort
Smooth, flared archslot	Improved sliding mechanics – less binding Less friction than traditional ceramics
Patented mushroom style base	Superior bond performance and strength Consistent debonding – greater confidence without sacrificing strength
Contoured base design	Optimum adaptation to tooth surfaces
Polycrystalline Ceramic	Unmatched clarity and improved fracture strength

Roth Prescription - Maxillary Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central	12°	5°	0°		R	Q07200	Q07300
	12°	5°	0°		L	Q07201	Q07301
Lateral	8°	9°	0°		R	Q07202	Q07302
	8°	9°	0°		L	Q07203	Q07303
Cuspid	0°	11°	4° MO		R	Q07204	Q07304
	0°	11°	4° MO		L	Q07205	Q07305
	0°	11°	4° MO	Distal	R	Q07206	Q07306
	2°	11°	4° MO	Distal	L	Q07207	Q07307
1st/2nd	-7°	0°	0°		R/L	Q07208	Q07308
Bicuspid	-7°	0°	0°	Center	R/L	Q07210	Q07310



Roth Prescription - Mandibular Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central/ Lateral	0°	0°	0°		R/L	Q07220	Q07320
Cuspid	-11°	5°	2° MO		R	Q07224	Q07324
	-11°	5°	2° MO		L	Q07225	Q07325
	-11°	5°	2° MO	Distal	R	Q07226	Q07326
	-11°	5°	2° MO	Distal	L	Q07227	Q07327
1st	-17°	0°	0°		R/L	Q07228	Q07328
Bicuspid	-17°	0°	0°	Center	R/L	Q07230	Q07330
2nd	-22°	0°	0°		R/L	Q07232	Q07332
Bicuspid	-22°	0°	0°	Center	R/L	Q07234	Q07334



Roth Prescription - Kits

Description	.018	.022
Roth FLI® Clear without hook, 5-5 U/L	KQ07200	KQ07300
Roth FLI® Clear hook on cuspid, 5-5 U/L	KQ07201	KQ07301
Roth FLI® Clear hook on cuspid/bicuspid, 5-5 U/L	KQ07202	KQ07302
Roth FLI® Clear hook on cuspid, 3-3 U	KQ07203	KQ07303
Roth FLI® Clear without hook, 3-3 U/L	KQ07204	KQ07304
Roth FLI® Clear hook on cuspid, 3-3 U/L	KQ07205	KQ07305

bracket systems

FLI® CLEAR BRACKETS: MBT PRESCRIPTION



MBT Prescription - Maxillary Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central	17°	4°	0°		R	Q07400	Q07500
	17°	4°	0°		L	Q07401	Q07501
Lateral	10°	8°	0°		R	Q07402	Q07502
	10°	8°	0°		L	Q07403	Q07503
Cuspid	-7°	8°	0°		R	Q07414	Q07514
	-7°	8°	0°		L	Q07415	Q07515
	-7°	8°	0°	Distal	R	Q07416	Q07516
	-7°	8°	0°	Distal	L	Q07417	Q07517
Extraction case	0°	8°	0°		R	Q07404	Q07504
	0°	8°	0°		L	Q07405	Q07505
	0°	8°	0°	Distal	R	Q07406	Q07506
	0°	8°	0°	Distal	L	Q07407	Q07507
1st/2nd	-7°	0°	0°		R/L	Q07408	Q07508
Bicuspid	-7°	0°	0°	Center	R/L	Q07410	Q07510



MBT Prescription - Mandibular Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central/ Lateral	-6°	0°	0°		R/L	Q07422	Q07522
Cuspid	-6°	3°	0°		R	Q07436	Q07536
	-6°	3°	0°		L	Q07437	Q07537
	-6°	3°	0°	Distal	R	Q07438	Q07538
	-6°	3°	0°	Distal	L	Q07439	Q07539
Extraction case	0°	3°	0°		R	Q07424	Q07524
	0°	3°	0°		L	Q07425	Q07525
	0°	3°	0°	Distal	R	Q07426	Q07526
	0°	3°	0°	Distal	L	Q07427	Q07527
1st	-12°	0°	0°		R/L	Q07428	Q07528
Bicuspid	-12°	0°	0°	Center	R/L	Q07430	Q07530
2nd	-17°	0°	0°		R/L	Q07432	Q07532
Bicuspid	-17°	0°	0°	Center	R/L	Q07434	Q07534

MBT Prescription - Kits

Description	.018	.022
MBT FLI® Clear without hook, 5-5 U/L	KQ07400	KQ07500
MBT FLI® Clear hook cuspid, 5-5 U/L	KQ07401	KQ07501
MBT FLI® Clear hook cuspid/bicuspid, 5-5 U/L	KQ07402	KQ07502
MBT FLI® Clear hook cuspid, 3-3 U	KQ07403	KQ07503
MBT FLI® Clear without hook, 3-3 U/L	KQ07404	KQ07504
MBT FLI® Clear hook cuspid, 3-3 U/L	KQ07405	KQ07505
Extraction Case	.018	.022
MBT FLI® Clear 0° T cuspid without hook, 5-5 U/L	KQ07406	KQ07506
MBT FLI® Clear 0° T cuspid hook cuspid, 5-5 U/L	KQ07407	KQ07507
MBT FLI® Clear 0° T cuspid hook cuspid/bicuspid, 5-5 U/L	KQ07408	KQ07508
MBT FLI® Clear 0° T cuspid hook cuspid, 3-3 U	KQ07409	KQ07509
MBT FLI® Clear 0° T cuspid without hook, 3-3 U/L	KQ07410	KQ07510
MBT FLI® Clear 0° T cuspid hook cuspid, 3-3 U/L	KQ07411	KQ07511

Ricketts® Original Biopressive® Prescription - Maxillary Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018
Central	22°	0°	0°		R/L	Q07600
Lateral	14°	8°	0°		R	Q07602
	14°	8°	0°		L	Q07603
Cuspid	7°	5°	0°		R	Q07604
	7°	5°	0°		L	Q07605
	7°	5°	0°	Distal	R	Q07606
	7°	5°	0°	Distal	L	Q07607
1st/2nd Bicuspid	0°	0°	0°	Center	R/L	Q07610



Ricketts® Original Biopressive® Prescription - Mandibular Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018
Central/ Lateral	0°	0°	0°		R/L	Q07620
Cuspid	7°	5°	0°		R	Q07624
	7°	5°	0°		L	Q07625
	7°	5°	0°	Distal	R	Q07626
	7°	5°	0°	Distal	L	Q07627
1st/2nd Bicuspid	0°	0°	0°	Center	R/L	Q07610



Ricketts® Original Biopressive® Prescription - Kits

Description	.018
Ricketts® Biopressive® FLI® Clear hook on cuspid/bicuspid, 5-5 U/L	KQ07602
Ricketts® Biopressive® FLI® Clear hook on cuspid, 3-3 U	KQ07603