


Altitude SL® Brackets

RMO®'s Altitude SL® Brackets are the new solution to low-friction, passive self-ligation technology that delivers trusted results and treatment efficiency. Fast and efficient treatment is attained by effortless door activation and optimal sliding mechanics. Altitude SL® offers reliable strength in construction combined with a smooth rounded profile for enhanced patient comfort.



MATERIAL: 17-4PH Stainless Steel Pad
BASE: 304 Stainless Steel 80-gauge Mesh
PRESCRIPTIONS: Roth & MBT™
SLOT SIZE: .018 & .022
TYPE: Passive Self-Ligation

FEATURES

Smooth, Rounded Profile
 Ball End Hook
 Optional Ligation
 Torque-in-Base
 Two-Piece Mesh Construction
 Smooth Rounded Slot Lead-in
 Patented Rail Design
 Door Notch
 Door Closed Gingivally
 Base Palmer Notation
 17-4PH Grade SS

BENEFITS

Enhanced patient comfort
 Added patient comfort through minimized cheek irritation
 Ideal finishing and ligation
 Optimal slot alignment reduces guesswork and chair time
 Consistent and predictable bonding
 Easy wire placement and insertion for quicker wire changes
 Lower friction and binding on wire
 Provides reliable, stable door retention
 Quick and easy to open and close
 No instrument required
 Keeps door closed positively with chewing
 Easy, permanent bracket identification
 Molded for high strength and rigidity

bracket systems

ALTITUDE SL® BRACKETS: ROTH PRESCRIPTION

Roth Prescription - Maxillary Brackets



Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central	12°	5°	0°		R		M11400
	12°	5°	0°		L		M11401
Lateral	8°	9°	0°		R		M11402
	8°	9°	0°		L		M11403
Cuspid	-2°	10°	0°	Distal	R		M11406
	-2°	10°	0°	Distal	L		M11407
1st/2nd	-7°	0°	0°	Distal	R		M11410
Bicuspid	-7°	0°	0°	Distal	L		M11411

Roth Prescription - Mandibular Brackets



Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central/ Lateral	-1°	0°	0°		R/L		M11420
Cuspid	-11°	7°	0°	Distal	R		M11426
	-11°	7°	0°	Distal	L		M11427
1st	-17°	0°	0°	Distal	R		M11430
Bicuspid	-17°	0°	0°	Distal	L		M11431
2nd	-22°	0°	0°	Distal	R		M11434
Bicuspid	-22°	0°	0°	Distal	L		M11435

Roth Prescription - Kits

Description	.018	.022
Roth Altitude SL®, with hook cuspid/bicuspid, 5-5		KM11402

MBT Prescription - Maxillary Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central	17°	4°	0°		R	M11000	M11100
	17°	4°	0°		L	M11001	M11101
Lateral	10°	8°	0°		R	M11002	M11102
	10°	8°	0°		L	M11003	M11103
Cuspid	-7°	8°	0°	Distal	R	M11006	M11106
	-7°	8°	0°	Distal	L	M11007	M11107
1st/2nd	-7°	0°	0°	Distal	R	M11010	M11110
Bicuspid	-7°	0°	0°	Distal	L	M11011	M11111



MBT Prescription - Mandibular Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central/ Lateral	-6°	0°	0°		R/L	M11020	M11120
Cuspid	-6°	3°	0°	Distal	R	M11026	M11126
	-6°	3°	0°	Distal	L	M11027	M11127
1st	-12°	2°	0°	Distal	R	M11030	M11130
Bicuspid	-12°	2°	0°	Distal	L	M11031	M11131
2nd	-17°	2°	0°	Distal	R	M11034	M11134
Bicuspid	-17°	2°	0°	Distal	L	M11035	M11135



MBT Prescription - Kits

Description	.018	.022
MBT Altitude SL®, with hook cuspid/bicuspid, 5-5	KM11002	KM11102