

FLI® Wire Orthonol® Nickel-Titanium Preformed Natural Aesthetic Coated Arch Wires

- Natural appearance that will last until your next wire change
- Maintains expected wire properties and performance characteristics of the wire
- Coating on labial surface only
- Nickel-titanium alloys are shape memory alloys
- Orthonol® is superelastic nickel-titanium in preformed natural arch shapes
- Provides an ideal continuous gentle force at a constant, uninterrupted force level for tooth movement
- Reduces chair time by extending time between visits

Shape	Size (inch)	Size (mm)	Package Contains	Maxillary	Mandibular
Round	0.012	(0.305)	10	A07011-A	A07031-A
	0.014	(0.356)	10	A07000-A	A07020-A
	0.016	(0.406)	10	A07001-A	A07021-A
	0.018	(0.457)	10	A07002-A	A07022-A
Square	.016 x .016	(0.406 x 0.406)	10	A07004-A	A07024-A
	.020 x .020	(0.508 x 0.508)	10	A07013-A	A07033-A
Rectangular	.016 x .022	(0.406 x 0.559)	10	A07005-A	A07025-A
	.017 x .025	(0.432 x 0.635)	10	A07006-A	A07026-A
	.018 x .025	(0.457 x 0.635)	10	A07008-A	A07028-A
	.019 x .025	(0.483 x 0.635)	10	A07009-A	A07029-A

FLI® Wire T.M.A. Preformed Natural Aesthetic Coated Arch Wires

- Molybdenum based alloy - nickel free

Shape	Size (inch)	Size (mm)	Package Contains	Maxillary	Mandibular
Rectangular	.016 x .022	(0.406 x 0.559)	10	A07125-A	A07134-A
	.017 x .025	(0.432 x 0.635)	10	A07126-A	A07135-A
	.019 x .025	(0.483 x 0.635)	10	A07128-A	A07137-A