









Micro Rapid Palatal Expander


- Uniform expansion with substantially smaller body than competing expanders.
- Features longer arms providing additional anterior support and stability to the patient
- The Micro Expander can also be used for mandible expansion
-  Cross sectional width 4.95mm
-  Cross sectional height 4.05mm
-  Arms diameter 1.5mm

Description		Maximum Expansion	Package Contains	Order Number
12mm	8mm	35 turns	1	A00670
16mm	11mm	50 turns	1	A00671






Orthogonal Rapid Palatal Expander


- Ideal for smaller or narrow palate applications
- Orthogonal arms allow fabrication of devices with limited dimensions while maintaining comfort in size and stability
-  Cross sectional width 5mm
-  Cross sectional height 4.5mm
-  Arms diameter 1.5mm

Description		Maximum Expansion	Package Contains	Order Number
12mm	8mm	40 turns	1	A00672
14mm	10mm	50 turns	1	A00673
16mm	12mm	60 turns	1	A00674



Pre-Bent Rapid Palatal Expander

- Longstanding reliability and performance for any lab application
- Arms pre-bent to precise anatomical angulation for ease of adaptation and fabrication
-  Cross sectional width 11mm
-  Cross sectional height 4mm
-  Arms diameter 1.5mm

Description		Maximum Expansion	Package Contains	Order Number
12mm	8mm	35 turns	1	A00676
14mm	9mm	40 turns	1	A00677
16mm	11mm	50 turns	1	A00678
18mm	13mm	60 turns	1	A00679