

### Facebow, Ricketts®

- Flexible, permits natural expansion
- Outer bow loop provides vertical auxiliary control
- Laser welded and brazed into a strong single unit



Description	Package Contains	Order Number
Non-extraction/Large Arch	1	J00400
Extraction/Small Arch	1	J00401

### Facebow, Penta-Form™ Loop-Style

- 5 sizes of inner bow; length in millimeters is measured from distal end of loop along inner arch to distal end of opposite loop
- .071" made entirely of biomedical spring stainless steel wire outer bow
- .045" inner bow
- Laser welded joint for enhanced strength and reliability
- Marked with size and midline indication



Size	Package Contains	Order Number
1 83mm	1	J01290
2 90mm	1	J01291
3 97mm	1	J01292
4 104mm	1	J01293
5 111mm	1	J01294
Universal Full Form	1	J01297

### Spring Gear®

- Comfortable, washable neck pads
- Continuously exerts the same force on the teeth regardless of head movement
- NiTi Coil spring to insure same force regardless of head movement



Description	Package Contains	Order Number
Spring-Gear®	1 Spring-Gear	J00521
Replacement Neck Pads	10	J00525

### Safety Headgear System

- Two-part safety release modules safely disengage when any excessive force is put on the face bow
- Easy for patient to reassemble after release
- Provides consistent traction force



Cervical Neckpad with 2 Modules	Package Contains	Order Number
Medium (14 oz.) white straps	1	J00390
Heavy (18 oz.) gray straps	1	J00391

High-Pull Head Strap with 2 Modules	Package Contains	Order Number
Medium (14 oz.) white straps	1	J00392
Heavy (18 oz.) gray straps	1	J00393

Description	Package Contains	Order Number
Neckpad	1	J00394
High-Pull Head strap	1	J00395
Modules Medium (14 oz.) white straps	2	J00396
Modules Heavy (18 oz.) gray straps	2	J00397