

intraoral & extraoral

ELASTOMERICS: ENERGY PAK™ ELASTICS



1/8" 3/16" 1/4" 5/16" 3/8" 1/2" 5/8" 3/4"



0.32cm 0.48cm 0.64cm 0.79cm 0.95cm 1.27cm 1.59cm 1.90cm

Force Values

General Purpose

Cuspid to Cuspid

Light	57gm (2oz.)	parachute J-01100	hang glider J-01101	balloon J-01102	dirigible J-01103	helicopter J-01104		skating J-01156	snow shoeing J-01157
Medium	85gm (3oz.)		motor home J-01111	truck J-01112	train J-01113	monorail J-01114			
Heavy	99gm (3.5oz.)	surfboard J-01120	canoe J-01121	sailboat J-01122	speedboat J-01123	submarine J-01124		toboggan J-01166	snow skiing J-01167
Heavy Plus	128gm (4.5oz.)		motorcycle J-01131	racecar J-01132	car J-01133	snowmobile J-01134			
Extra Heavy	142gm (5oz.)	jumprope J-01140	skateboard J-01141	unicycle J-01142	bicycle J-01143	jogging J-01144			
Assorted Neon									
Neon	99gm (3.5oz.)	convertible J-01400	roller skate J-01401	moped J-01402	jet ski J-01403				
Neon	128gm (4.5oz.)	taxi J-01404	snowboard J-01405	van J-01406	walking J-01407				
Neon	170gm (6.0oz)	scooter J-01408	bobsled J-01409	police J-01410	paddle board J-01411				
Extraoral									
Light	454gm (16oz.)		bus J-01171	ocean liner J-01172	jumbo jet J-01173	spaceship J-01174	rocket J-01175		
Heavy	908gm (32oz.)		rickshaw J-01181	horse J-01182	chariot J-01183	covered wagon J-01184	stage coach J-01185		

- J01100-J01144: 100 elastics/bag; 50 bags/box
- J01400-J01411: 100 elastics/bag; 50 bags/box
- J01156-J01167: 50 elastics/bag; 40 bags/box
- J01171-J01185: 50 elastics/bag; 20 bags/box

- **LF** Available in latex free
- Once pull equals three times slack diameter
- To achieve ounce pull, two elastics must be used (only applies to extraoral elastics)