

1st Molar - Maxillary Single Tubes



Description	Style	Hook	R/L	.018	.022
Standard Edgewise, without distal extension	Weld		R/L		A00210
Standard Edgewise, without distal extension	Bond		R/L		A08012



Standard Edgewise, distal extension	Weld		R/L	A00160	A00101
Standard Edgewise, distal extension	Bond		R/L	A08000	A08010



Standard Edgewise, distal extension long flange	Weld		R/L	A00952	A00953
Standard Edgewise, distal extension long flange	Weld	Mesial	UR/LL	A00954	A00955
Standard Edgewise, distal extension long flange	Weld	Mesial	UL/LR	A00956	A00957
Standard Edgewise, distal extension long flange	Bond	Mesial	UR/LL	A08006	A08016
Standard Edgewise, distal extension long flange	Bond	Mesial	UL/LR	A08007	A08017



FLI® Standard Edgewise	Weld	Mesial	R	A06770	A06774
FLI® Standard Edgewise	Weld	Mesial	L	A06771	A06775
FLI® Standard Edgewise	Bond	Mesial	R	A08770	A08774
FLI® Standard Edgewise	Bond	Mesial	L	A08771	A08775



Long flange, distal offset	Weld		R/L		A00975
Long flange, distal offset	Weld	Mesial	UR/LL		A00977
Long flange, distal offset	Weld	Mesial	UL/LR		A00979



FLI® Tube	Weld	Mesial	R	A06730	A06734
FLI® Tube	Weld	Mesial	L	A06731	A06735
FLI® Tube	Bond	Mesial	R	A08730	A08734
FLI® Tube	Bond	Mesial	L	A08731	A08735







Andrews, convertible	Weld	Mesial	R	A06112	A06152
Andrews, convertible	Weld	Mesial	L	A06113	A06153
Andrews, convertible	Bond	Mesial	R	A08412	A08452
Andrews, convertible	Bond	Mesial	L	A08413	A08453



FLI® One Piece	Bond	Mesial	R	A18750	A18754
FLI® One Piece	Bond	Mesial	L	A18751	A18755

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	6° D							
0°	0°	6° D							
0°	0°	6° D							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	7° D							
-10°	0°	7° D							
-10°	0°	7° D							
-10°	0°	7° D							
-10°	0°	10° D							
-10°	0°	10° D							

1st Molar - Maxillary Single Tubes

Description	Style	Hook	R/L	.018	.022
 Bench Bioprogressive®, convertible Bench Bioprogressive®, convertible Bench Bioprogressive®, convertible Bench Bioprogressive®, convertible	Weld	Mesial	R	A01194	A01294
	Weld	Mesial	L	A01195	A01295
	Bond	Mesial	R	A08294	A08494
	Bond	Mesial	L	A08295	A08495
 Roth, convertible Roth, convertible Roth, convertible Roth, convertible	Weld	Mesial	R	A05664	A05692
	Weld	Mesial	L	A05665	A05693
	Bond	Mesial	R	A08634	A08636
	Bond	Mesial	L	A08635	A08637
 FLI® MBT FLI® MBT FLI® MBT FLI® MBT	Weld	Mesial	R	A05422	A05522
	Weld	Mesial	L	A05423	A05523
	Bond	Mesial	R	A08562	A08566
	Bond	Mesial	L	A08563	A08567
 FLI® Roth FLI® Roth FLI® Roth FLI® Roth	Weld	Mesial	R	A05600	A05602
	Weld	Mesial	L	A05601	A05603
	Bond	Mesial	R	A08200	A08202
	Bond	Mesial	L	A08201	A08203

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							