


1st Molar - Maxillary Double Tubes

Description	Style	Hook	R/L	.018	.022
 FLI® Roth	Weld	Mesial	R	A05616	A05618
FLI® Roth	Weld	Mesial	L	A05617	A05619
FLI® Roth	Bond	Mesial	R	A08216	A08218
FLI® Roth	Bond	Mesial	L	A08217	A08219
Roth, convertible	Weld	Mesial	R	A05650	A05670
Roth, convertible	Weld	Mesial	L	A05651	A05671
Roth, convertible	Bond	Mesial	R	A08600	A08620
Roth, convertible	Bond	Mesial	L	A08601	A08621
Roth, convertible	Weld	Mesial	R	A05652	A05672
Roth, convertible	Weld	Mesial	L	A05653	A05673
Roth, convertible	Bond	Mesial	R	A08602	A08622
Roth, convertible	Bond	Mesial	L	A08603	A08623
Roth, convertible	Weld	Mesial	R	A05684	A05686
Roth, convertible	Weld	Mesial	L	A05685	A05687
Roth, convertible	Bond	Mesial	R	A08684	A08686
Roth, convertible	Bond	Mesial	L	A08685	A08687
Roth, without cap	Weld	Mesial	R	A05688	A05690
Roth, without cap	Weld	Mesial	L	A05689	A05691
Roth, without cap	Bond	Mesial	R	A08688	A08690
Roth, without cap	Bond	Mesial	L	A08689	A08691

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-14°	0°	-14° D	-14°	0°	14° D				
-14°	0°	-14° D	-14°	0°	14° D				
-14°	0°	-14° D	-14°	0°	14° D				
-14°	0°	-14° D	-14°	0°	14° D				
-14°	0°	14° D				0° D	14° D	Gingival	.045
-14°	0°	14° D				0° D	14° D	Gingival	.045
-14°	0°	14° D				0° D	14° D	Gingival	.045
-14°	0°	14° D				0° D	14° D	Gingival	.045
-14°	0°	14° D				0° D	14° D	Occlusal	.045
-14°	0°	14° D				0° D	14° D	Occlusal	.045
-14°	0°	14° D				0° D	14° D	Occlusal	.045
-14°	0°	14° D				0° D	14° D	Occlusal	.045
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				