

bands

TRUFORM BANDS MAXILLARY FIRST MOLAR RIGHT AND LEFT

- Anatomical design – right and left shape for precise and comfortable fitting
- Size identification and Palmer notation are permanently marked to withstand heat sterilization
- Accurate sizes – proportionately graduated sizes of superior RMO temper for use on primary molars and first permanent molars
- Dimensional strength remains stable without deformation during trial fitting
- Withstands mastication and external appliance stress



Primary Second Molar and First Permanent Molar, Equal Number of Rights and Lefts

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total	
Kit Contains	4	4	4	4	4	4	4	4	4	4	4	4	6	6	8	10	12	12	12	10	8	8	8	6	6	6	4	4	4	4	4	4	4	4	4	4	4	214
Introductory Kit Order Number																																					K00518	



First Permanent Only, Equal Number of Rights and Lefts

Band Size	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total										
Kit Contains	4	6	6	8	8	10	12	12	12	10	10	8	8	6	6	6	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	180
Introductory Kit Order Number																												K00520									
Band Size	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	Total																		
Kit Contains	2	2	2	4	4	4	6	8	10	10	10	10	8	6	6	4	2	2	100																		
Mini Kit Order Number																		K00158																			

- Anatomical design – right and left shape for precise and comfortable fitting
- Size identification and Palmer notation are permanently marked to withstand heat sterilization
- Anatomically shaped – design conforms to buccal/lingual inclination and adapts accurately to the morphology of the tooth anatomy
- Accurate sizes – proportionately graduated sizes of superior RMO temper for use on primary molars and first permanent molars
- Dimensional strength remains stable without deformation during trial fitting
- Withstands mastication and external appliance stress



Primary Second Molar and First Permanent Molar, Equal Number of Rights and Lefts

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total	
Kit Contains	4	4	4	4	4	4	4	4	4	4	6	8	8	10	12	12	12	12	10	8	6	6	6	6	4	4	4	4	4	4	4	4	4	4	4	4	4	218
Introductory Kit Order Number																																					K00558	



First Permanent Only, Equal Number of Rights and Lefts

Band Size	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total												
Kit Contains											8	8	10	12	12	12	12	12	8	8	6	6	6	4	4	4	4	4	4	4	4	4	4	4	4	4	4	176
Introductory Kit Order Number																									K00560													
Band Size	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	Total																			
Kit Contains	2	2	2	4	4	4	6	8	10	10	10	10	8	6	6	4	2	2	100																			
Mini Kit Order Number																			K00159																			

bands

TRUFORM BANDS MAXILLARY FIRST MOLAR UNIVERSAL

- Universally shaped – functional fit with complete adaptation to tooth anatomy
- Size identification and Palmer notation are permanently marked to withstand heat sterilization
- Accurate sizes – proportionately graduated sizes of superior RMO temper for use on primary molars and first permanent molars
- Buccal groove for orientation to tooth anatomy
- Dimensional strength remains stable without deformation during trial fitting
- Withstands mastication and external appliance stress



Primary Second Molar and First Permanent Molar

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total				
Band Size																																									
Kit Contains	2	2	2	2	2	2	2	3	3	3	3	4	4	5	7	8	9	9	9	9	8	7	6	6	5	5	4	4	3	3	2	2	2	2	2	2	2	2	153		
Introductory Kit Order Number																																				K00510					
Mini Kit Contains	2	2	2	2	2	2	2	2	2	2	2	2	3	3	4	5	6	6	6	5	4	4	4	3	3	3	2	2	2	2	2	2	2	2	2	2	2	2	107		
Mini Kit Order Number																																				K00050					

Primary Second Molar and First Permanent Molar, packaged for Equal Number of Rights and Lefts for Tube/Bracket Placement

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total				
Band Size																																									
Kit Contains	4	4	4	4	4	4	4	4	4	4	4	4	6	6	8	10	12	12	12	12	10	8	8	8	6	6	6	4	4	4	4	4	4	4	4	4	4	4	214		
Introductory Kit Order Number																																				K00512					

First Permanent Only

	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total
Band Size																										
Kit Contains	4	5	6	7	8	10	10	10	9	8	7	6	6	6	5	4	4	4	3	2	2	2	2	2	134	
Introductory Kit Order Number																								K00514		
Band Size	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28								
Kit Contains	2	2	2	4	4	4	6	7	10	10	10	8	6	6	4	2	2								100	
Mini Kit Order Number																		K00148								

First Permanent Only, packaged for Equal Number of Rights and Lefts for Tube/Bracket Placement

	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total
Band Size																										
Kit Contains	6	6	8	8	10	12	12	12	12	10	10	8	8	8	6	6	6	4	4	4	4	4	4	4	180	
Introductory Kit Order Number																								K00516		

- Universally shaped – functional fit with complete adaptation to tooth anatomy
- Size identification and Palmer notation are permanently marked to withstand heat sterilization
- Accurate sizes – proportionately graduated sizes of superior RMO temper for use on primary molars and first permanent molars
- Lingual groove for orientation to tooth anatomy
- Developed through years of RMO band manufacturing
- Dimensional strength remains stable without deformation during trial fitting
- Withstands mastication and external appliance stress



Primary Second Molar and First Permanent Molar

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total		
Kit Contains	2	2	2	2	2	2	2	2	3	4	5	6	7	8	9	9	9	9	8	7	5	5	4	4	3	3	2	2	2	2	2	2	2	2	2	2	2	2	151
Introductory Kit Order Number																																				K00550			
Mini Kit Contains	2	2	2	2	2	2	2	2	2	3	4	4	5	6	6	6	6	5	4	3	3	3	3	2	2	2	2	2	2	2	2	2	2	2	2	2	2	109	
Mini Kit Order Number																																				K00051			

Primary Second Molar and First Permanent Molar, packaged for Equal Number of Rights and Lefts for Tube/Bracket Placement

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total	
Kit Contains	4	4	4	4	4	4	4	4	4	6	8	8	10	12	12	12	12	12	10	8	6	6	6	6	4	4	4	4	4	4	4	4	4	4	4	4	4	218
Introductory Kit Order Number																																				K00552		

First Permanent Only

Band Size												12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total														
Kit Contains												6	7	8	10	10	10	10	9	7	6	5	4	4	4	3	3	2	2	2	2	2	2	2	2	2	132														
Introductory Kit Order Number																																				K00554															
Band Size												11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28																		Total				
Kit Contains												2	2	2	4	4	4	6	8	10	10	10	10	8	6	6	4	2	2																						100
Mini Kit Order Number																																				K00149															

First Permanent Only, packaged for Equal Number of Rights and Lefts for Tube/Bracket Placement

Band Size												12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	Total		
Kit Contains												8	8	10	12	12	12	12	12	8	8	6	6	6	4	4	4	4	4	4	4	4	4	4	4	4	4	4	176
Introductory Kit Order Number																																				K00556			

bands

TRUFORM BANDS MAXILLARY/MANDIBULAR SECOND MOLAR RIGHT AND LEFT

- Anatomical right and left shape with pre-positioned buccal groove provides positive retention with minimum effort and adaptation
- Size identification and Palmer notation are permanently marked to withstand heat sterilization
- Anatomically shaped – design conforms to buccal/lingual inclination and adapts to morphology of the tooth anatomy
- Dimensional strength remains stable without deformation during trial fitting
- Withstands mastication and external appliance stress
- Full complement of sizes – 26 size increments of superior RMO temper



Maxillary Second Molar, packaged for Equal Number of Rights and Lefts

Band Size	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Total
Kit Contains	4	4	4	4	4	6	6	8	8	10	12	12	12	12	10	8	8	8	8	6	6	4	4	4	4	4	180

Introductory Kit Order Number **K00522**



Mandibular Second Molar, packaged for Equal Number of Rights and Lefts

Band Size	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Total
Kit Contains	4	4	4	4	6	8	8	10	12	12	12	12	12	10	8	6	6	6	6	4	4	4	4	4	4	4	178

Introductory Kit Order Number **K00564**



- Universally shaped – functional fit with complete adaptation to tooth anatomy; gingival festoon eliminates tissue impingement
- Full complement of sizes – choice of correct size for maximum retention at desired position
- Size identification – color-coded black with letter for size selections and maxillary identification
- RMO-engineered ductility and precise form adapts to the tooth contours with secure retention



Universal Bracket Placement

Band Size	G	G+	H	H+	I	I+	J	J+	K	K+	L	L+	M	M+	N	N+	Total
Kit Contains	2	3	4	5	7	10	13	16	17	15	13	10	7	4	3	2	131
Introductory Kit Order Number	K00329																

Packaged for Equal Number of Rights and Lefts for Bracket Placement

Band Size	G	G+	H	H+	I	I+	J	J+	K	K+	L	L+	M	M+	N	N+	Total
Kit Contains	4	4	6	8	10	12	16	20	22	18	16	12	8	6	4	4	170
Introductory Kit Order Number	K00330																

bands

TRUFORM BANDS MANDIBULAR FIRST AND SECOND BICUSPID



- Universally shaped – functional fit with complete adaptation to tooth anatomy
- Correct band width – secure anchorage is provided without gingival extension
- Full complement of sizes – choice of correct size for maximum retention at desired position
- Size identification – color-coded red with letter for size selections and mandibular identification
- RMO-engineered ductility and precise form adapts to the tooth contours with secure retention; festooned gingival

Mandibular First Bicuspid, Universal Bracket Placement

Band Size	F	F+	G	G+	H	H+	I	I+	J	J+	K	K+	L	L+	M	M+	Total
Kit Contains	2	3	6	8	9	11	12	13	10	7	5	3	2	2	2	2	97
Introductory Kit Order Number	K00331																

Mandibular First Bicuspid, packaged for Equal Number of Rights and Lefts for Bracket Placement

Band Size	F	F+	G	G+	H	H+	I	I+	J	J+	K	K+	L	L+	M	M+	Total
Kit Contains	4	4	8	10	12	14	16	16	14	10	6	4	4	4	4	4	134
Introductory Kit Order Number	K00332																



- Apron contoured tightly to buccal surface and adapts smoothly with mesial and distal angles
- Size Identification - color-coded red with number for size selections and mandibular identification

Mandibular Second Bicuspid, Universal Bracket Placement

Band Size	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	Total	
Kit Contains	2	2	2	3	4	4	4	4	5	6	7	8	9	9	8	7	6	5	4	3	2	2	2	2	2	2	2	123	
Introductory Kit Order Number	K00333																												

Mandibular Second Bicuspid, packaged for Equal Number of Rights and Lefts for Bracket Placement

Band Size	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	Total	
Kit Contains	4	4	4	4	6	6	6	6	6	8	8	10	12	12	12	10	10	8	6	6	4	4	4	4	4	4	4	4	176
Introductory Kit Order Number	K00334																												



- Universally shaped – apron bands permit offset incisal placement of bracket
- Apron contours tightly to labial surface and adapts smoothly with mesial and distal angles
- Rolled gingival and incisal edges provide total retention and minimize cement ridge build-up
- Specific incisal-gingival taper for maxillary canines with added ductility for superior fitting qualities
- Full complement of sizes – clinically tested for maximum utilization
- Size identification – color-coded black with letter for size selections and maxillary identification
- RMO-engineered ductility and precise form adapts to the tooth contours with secure retention

Maxillary Apron Cuspid, Universal Bracket Placement

Band Size	D	D+	E	E+	F	F+	G	G+	H	H+	I	I+	J	J+	Total
Kit Contains	2	3	4	6	9	11	13	13	11	8	5	4	3	3	95
Introductory Kit Order Number	K00313														

Maxillary Apron Cuspid, packaged for Equal Number of Rights and Lefts for Bracket Placement

Band Size	D	D+	E	E+	F	F+	G	G+	H	H+	I	I+	J	J+	Total
Kit Contains	4	4	6	8	12	14	16	16	14	10	8	6	4	4	126
Introductory Kit Order Number	K00314														



- Size identification – color-coded red with letter for size selections and mandibular identification

Mandibular Apron Cuspid, Universal Bracket Placement

Band Size	B	B+	C	C+	D	D+	E	E+	F	F+	G	G+	H	H+	Total
Kit Contains	2	2	3	7	11	12	14	14	11	8	5	2	2	2	95
Introductory Kit Order Number	K00317														

Mandibular Apron Cuspid, packaged for Equal Number of Rights and Lefts for Bracket Placement

Band Size	B	B+	C	C+	D	D+	E	E+	F	F+	G	G+	H	H+	Total
Kit Contains	4	4	4	8	14	16	18	18	14	10	6	4	4	4	128
Introductory Kit Order Number	K00318														