

How to Seat RMO® Bands

There are two instruments that should be used when seating RMO® bands. They are Band Adapter (i00064) and the Band Pusher (i00300).

Band Adapter (RMO® i00064)



Durable plastic bite block and RMO® tempered stainless band pin

For bite or hand rocking bands

White plastic handle

Serrated pin

Replacement tip (RMO® i00065) (can be installed chairside)

Band Pusher (RMO® i00300)



Serrated stainless steel tip, non-removable

Stainless steel handle - square with chamfered corners

Tip end has square serrations

