

Make the right impression every time

with RMO®'s

Aluminum Impression Trays

RMO®'s Aluminum Impression Trays deliver flexible adjustability that no steel or plastic tray can match. With quality construction and colorful matte finishes, RMO®'s Aluminum Impression Trays adapt to every patient every time.

- Bendable and pliable for accurate impressions of any arch size
- Sterilizable
- 5 Maxillary/Mandibular sizes
- Color coded for rapid size identification
- Smooth rounded edges for optimal patient comfort
- Perforated design enables excellent alginate retention



MAXILLARY – package of 1		(width x length)	
J00304	Size 1	Pink	57 x 48 mm
J00305	Size 2	Blue	62 x 54 mm
J00306	Size 3	Red	68 x 60 mm
J00307	Size 4	Gold	75 x 65 mm
J00308	Size 5	Violet	84 x 70 mm

MANDIBULAR – package of 1			
J00364	Size 1	Pink	62 x 44 mm
J00365	Size 2	Blue	70 x 48 mm
J00366	Size 3	Red	75 x 50 mm
J00367	Size 4	Gold	82 x 53 mm
J00368	Size 5	Violet	87 x 55 mm

ACCESSORIES

- J00370 Stainless Steel Aluminum Impression Tray Tower
- J00371 Humidity/Storage Case (For Aluminum Impression Tray Tower)

KITS

- K00566 Maxillary and Mandibular Trays (1 of Each Size, 10 Total)
- K01564 Humidity/Storage Case & Tray Tower (No Trays Included)
- K01566 Humidity/Storage Case, Tray Tower, and 10 Trays (Maxillary and Mandibular - 1 of Each Size)

Stainless Steel Aluminum Impression Tray Tower

- Conveniently holds trays, impressions, or poured stone models
- Sterilizable, and reduces storage space in practice



Humidity/Storage Case (For Aluminum Impression Tray Tower)

- Store trays/impressions in a high humidity environment for less distortion before pouring models
- Humidity helps prevent crack formations
- Case is multi-functional and can be used for cold sterilization as well

