

intraoralextraoral

BULK ELASTICS, ENERGY CHAIN™

Bulk Elastics

- Finest quality latex surgical tubing (contains natural rubber latex which may cause allergic reactions)
- Top performance, fatigue life and force reliability
- Uniform precision cut for size and dimension
- May be used as separating elastics prior to band fitting


Inch	Size (mm)	Medium	Heavy	Package Contains
1/8"	(3.18)	J00129		approx. 2,000
3/16"	(4.76)	J00130		approx. 2,000
1/4"	(6.35)	J00131	J00134	approx. 2,000
5/16"	(7.94)	J00132	J00135	approx. 2,000



Plastic Pocket Containers

- Zippered polyethylene bag
- Elastic placement diagram on package
- Convenient pocket-size dispenser

Order Number	Package Contains
J00092	100




Energy Chain™ Elastics

- Latex-Free
- Tough and tasteless, non-toxic, non-staining, more consistent force delivery
- Easy-handling dispensing spools
- Spacing variation between modules allows proper selection, retraction, ligating, maxillary/mandibular/interarch application, etc.
- Available in clear, gray and an assortment of popular colors
- Package of one spool contains seven feet, six inches (2.3m)



	Clear	Gray	Blue	Green	Red	Purple	Package Contains
Closed	J00120	J00166	J00223	J00227	J00231	J00235	1 spool
Reduced	J00269	J00270	J00271	J00272	J00273	J00274	1 spool
Narrow	J00121	J00167	J00224	J00228	J00232	J00236	1 spool
Medium	J00122	J00168					1 spool




	Orange	Porsche Red	Black	Light Pink	Yellow	Package Contains
Closed	J00239	J00261	J00263	J00265	J00267	1 spool
Reduced	J00275	J00276	J00277	J00278	J00279	1 spool
Narrow	J00240	J00262	J00264	J00266	J00268	1 spool

Energy Chain™ Dispenser

- Organize the various colors and sizes of the Energy Chain spools
- Tinted plastic with a removable crossbar rod and sliding lid
- The Energy Chain is simply inserted through the 10 holes on the front portion of the dispenser box
- 3.25"H x 5"W x 3"D

Order Number	Package Contains
J00600	1 Dispenser



INTRAORALEXTRAORAL