

adhesives

RMO® SELF ETCHING LIGHT/DUAL CURE SYSTEMS/
RMO® BAND CEMENT / FLI® CERAMIC ADHESIVES

MonoLok™2 Self-Etching Light Cure Overview


- One-step etchant and bonding resin in one solution
- Time saving: eliminates one step in bonding procedure
- Convenient, easy to use dispensing system
- Economical and fast bonding/re-bonding of brackets
- Compatible with all light-cure adhesives



	Order Number	Package Contains
	K00211	1 Kit
Kit includes:		
0.1 ml single dose units		48
Micro brush tips		50
Brush holder		1

MonoLok™2 Self-Etching Dual Cure Overview

- Unique, easy-to-use dispensing system
- Two components of mixture are separated by foil that is easy to penetrate with a supplied microbrush
- Easy to mix
- One-step etchant and bonding resin in mixed solution to prepare tooth surface
- Time saving, to prepare tooth surface
- Dual cure: chemical reaction and/or visible light



	Order Number	Package Contains
	K00210	1 Kit
Kit includes:		
2 x 0.1 ml Single Mix units		35
Micro brushes		35

RMO® Band Cement

	Order Number	Package Contains
RMO® Band Cement 5gm Syringe	J03280	1



FLI® Ceramic Light Activated Adhesive

	Order Number	Package Contains
	K00200	1 Kit
Kit includes:		
FLI Ceramic Adhesive 3.5gm Syringe		
Sealant Resin 3 ml bottle		



ADHESIVES