

## Ligatures and Separators



- Convenient dispensing
- Strong, tasteless, maintains retention


	Gray	Clear	Color	Package Contains
<b>Ligature "S" Shaped Dispenser</b> .123" (3.12mm) outside diameter (10 dispensers)	<b>J00156</b>	<b>J00154</b>		1,000

### Synergy® Low-Friction Ligatures

- Silicone injected
- Reduced friction
- Good stain resistance
- Low profile

.123" (3.12mm) outside diameter  
(10 dispensers)

			<b>J00151</b> /white opaque	1,000
--	--	--	-----------------------------	-------



**Mini Single Case Ligature Stick**  
.121" (3.07mm) outside diameter  
100 sticks (10 ligatures per stick)

	J00350	J00345		1,000
			<ul style="list-style-type: none"> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #e0e0e0; border: 1px solid black; margin-right: 5px;"></span> <b>J00337</b>/pearl</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #a0c0ff; border: 1px solid black; margin-right: 5px;"></span> <b>J00338</b>/pearl blue</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #a0ffa0; border: 1px solid black; margin-right: 5px;"></span> <b>J00339</b>/pearl green</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #ffff00; border: 1px solid black; margin-right: 5px;"></span> <b>J00340</b>/fluorescent yellow</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #ff00ff; border: 1px solid black; margin-right: 5px;"></span> <b>J00341</b>/fluorescent pink</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #00ff00; border: 1px solid black; margin-right: 5px;"></span> <b>J00342</b>/fluorescent green</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #ff0000; border: 1px solid black; margin-right: 5px;"></span> <b>J00343</b>/fire engine red</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #0000ff; border: 1px solid black; margin-right: 5px;"></span> <b>J00344</b>/orchid</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #808080; border: 1px solid black; margin-right: 5px;"></span> <b>J00346</b>/smoke</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #add8e6; border: 1px solid black; margin-right: 5px;"></span> <b>J00347</b>/sky blue</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #ffa500; border: 1px solid black; margin-right: 5px;"></span> <b>J00348</b>/orange</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #000080; border: 1px solid black; margin-right: 5px;"></span> <b>J00349</b>/navy blue</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #0000ff; border: 1px solid black; margin-right: 5px;"></span> <b>J00351</b>/blue</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #008080; border: 1px solid black; margin-right: 5px;"></span> <b>J00352</b>/azure</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #006400; border: 1px solid black; margin-right: 5px;"></span> <b>J00353</b>/dark green</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #008000; border: 1px solid black; margin-right: 5px;"></span> <b>J00354</b>/green</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #ff0000; border: 1px solid black; margin-right: 5px;"></span> <b>J00355</b>/red</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #ffb6c1; border: 1px solid black; margin-right: 5px;"></span> <b>J00356</b>/light pink</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #dc143c; border: 1px solid black; margin-right: 5px;"></span> <b>J00357</b>/fuschia</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #d4af37; border: 1px solid black; margin-right: 5px;"></span> <b>J00358</b>/gold</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #808080; border: 1px solid black; margin-right: 5px;"></span> <b>J00359</b>/silver</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #e0e0e0; border: 1px solid black; margin-right: 5px;"></span> <b>J00360</b>/white</li> <li><span style="display: inline-block; width: 10px; height: 10px; background-color: #000000; border: 1px solid black; margin-right: 5px;"></span> <b>J00361</b>/black</li> </ul>	

### Mini Stick Colored Ligature Kit

Each kit contains 24 colors, with each compartment containing 1,000 ligatures. The lid contains a transparency which shows each part number, color, and position. Also included is a laminated color wheel which shows the 24 different colors allowing patients to choose their color chairside without touching the ligatures.



	Order Number	Package Contains
	<b>K00340</b>	24,000 ligatures
24 Compartment Box	<b>J00780</b>	1




# intraoralextraoral

LIGATURES/SEPARATORS, ROTATION WEDGES

## Ligatures and Separators

- Convenient dispensing
- Strong, tasteless, maintains retention

	Gray	Clear	Color	Package Contains
<b>Separator Sticks</b> .191" (4.85mm) outside diameter (100 sticks)			<b>J00177</b> /blue <b>J01451</b> /radio opaque	1,000
<b>Separators Cut</b> .170" (4.32mm) outside diameter			<b>J00178</b> /blue <b>J01450</b> /radio opaque	1,000 (approx.)
<b>Single Case Ligature Stick</b> .123" (3.12mm) outside diameter 45 sticks (22 ligatures per stick)	<b>J00221</b>	<b>J00220</b>		990
<b>Ligature Ringlets on Rings</b> .123" (3.12mm) outside diameter (10 rings)	<b>J00175</b>	<b>J00176</b>		1,000



## Elastomeric Ligatures

- Preloaded cane dispensers
- Easily engages under tie wings
- Retains resiliency over long periods of time
- Resists fluid absorption and odor



	Gray	Clear	Package Contains
<b>Ligature Canes</b> - .110" (2.79 mm) 20 canes of one color	<b>J00182</b>	<b>J00184</b>	1,000
<b>Ligature Canes</b> - .120" (3.05 mm) 20 canes of one color	<b>J00183</b>	<b>J00185</b>	1,000

## Elastomeric Rotation Wedges

- Produces rotational force by snapping perforated tab over tie wings of bracket, inserting arch wire in groove of wedge and ligating opposite tie wings of bracket to arch wire
- Wedge acts as fulcrum between bracket and arch wire
- Gray color



	Gray	Package Contains
	<b>J00174</b>	50 (5 sticks)

INTRAORALEXTRAORAL