




The ultimate alternative solution to standard model impression materials

TruForm™ Model Material is the ideal alternative solution to standard model impression materials by offering fast preparation and set, accurate patient arch forms, and convenient chairside efficiency.



APPLICATIONS AND USES:

- Indirect bonding (IDB) models ~ same day start appointments
- Same visit retainers ~ models can go home with patient
- Laboratory appliance fitting and installation
- Quick preparation consultation, diagnostic, and study models
- Better than traditional Stone or Plaster models
- Scan models with ease

FEATURES:

BENEFITS:

Quick set time (3-5 min)	No overnight drying – immediate appliance/retainer fabrication
	Quick chairside preparation – immediate use
Single use packet	Easy to mix and dispense
	Reduced waste
	No mess or mixing bowl required
Patented formula	No trimming or cleaning needed
	Smoother surface finish and consistency – reduced bubbles and porosity
Higher durability	Easy to maintain and store
	Less prone to breakage, chips, damage
	Reusable – no need to recreate model for repeated use

PART:

DESCRIPTION:

J03400	TruForm™ Model Material - 5 Single Use Packets (20 ml)
	1 arch per packet

P00776

