

## Tough O Energy Ligatures

- Trusted RMO® Energy Chain™ formula provides added strength and rebound characteristics
- Higher force retention than standard ligature ties keep wire in the slot longer
- Individual patient configuration – 20 ligatures per stick



<b>Material:</b>	Medical Grade Polyurethane (Latex Free)	
<b>Size:</b>	Thickness .032" (0.8 mm) Diameter .120" (3.0 mm)	
<b>Quantity:</b>	30 Sticks x 20 ligatures 600 ligatures per package	
<b>Color:</b>	Clear	
<b>Competition:</b>	Standard Elastomeric Ligatures Dentsply GAC Sani-Tie® American Orthodontics Unisticks Ortho Technology Power Sticks™	Ormco Molded "O"s 3M Unitek Mini-Stik™ Henry Schein Ortho O-Tie™

Features	Benefits
Proprietary formula * Same as RMO Energy Chain™	Tear and stain resistant
	Superior elasticity and reduced permanent set
Higher force retention	Increased torque expression of wire over extended time Wire stays in the slot longer
	No loss of friction - excellent retention during finishing phase
Increased elasticity	Increased interval time between appointments – fewer patient visits and freed up chair time
Medical grade polyurethane	Latex free – no allergy concerns
Grip tab design	Easy to hold and dispense
Individual patient configuration	20 ligatures – eliminates potential patient cross contamination

Part	Description
J00138	Tough O Energy Ligatures – 30 sticks x 20 ligatures

