

bands

RMO® TRUFORM BAND: MAXILLARY/MANDIBULAR SECOND MOLAR

- Anatomical right and left shape with pre-positioned buccal groove provides positive retention with minimum effort and adaptation
- Size identification and Palmer notation are permanently marked to withstand heat sterilization
- Anatomically shaped – design conforms to buccal/lingual inclination and adapts to morphology of the tooth anatomy
- Dimensional strength remains stable without deformation during trial fitting
- Withstands mastication and external appliance stress
- Full complement of sizes (26 sizes)

Second Maxillary Molar Right and Left



Band Size	Kit Contents
7	4
8	4
9	4
10	4
11	4
12	6
13	6
14	8
15	8
16	10
17	12
18	12
19	12
20	12
21	10
22	8
23	8
24	8
25	8
26	6
27	6
28	4
29	4
30	4
31	4
32	4
Total	180
Order#	K00522

Second Mandibular Molar Right and Left



Band Size	Kit Contents
7	4
8	4
9	4
10	4
11	6
12	8
13	8
14	10
15	12
16	12
17	12
18	12
19	12
20	10
21	8
22	6
23	6
24	6
25	6
26	4
27	4
28	4
29	4
30	4
31	4
32	4
Total	178
Order#	K00564