

bands

ACCESSORIES: PREWELDED AUXILIARIES

**Band Seating Lug**

- Used to seat preformed bands
- Furnished in 10-lug strips for easy welding to lingual

Description	Package Contains	Order Number
Bicuspid	100	A00718
Molar	100	A00719

**Lug/Cleat (Debnam)**

- Used as band seating lug and mesial-distal elastic hook
- Hook is bent out from band surface to engage elastics

Description	Package Contains	Order Number
Bicuspid	approx. 100	A00476
Molar	approx. 100	A00720

**Eyelets**

- Large, accepts wire and small diameter ligature thread

Description	Package Contains	Order Number
Regular	100	A00039

**Lingual Hook**

- Flat weldable base with smooth ball end for placement on molars

Description	Package Contains	Order Number
Regular	10	A00297

**Lingual Button**

- Exclusive round design allows for placement anywhere on band with minimum protrusion
- Can be used as a combination hook/band seating lug
- Ultra-smooth finish for improved patient comfort
- Easily welded to stainless steel bands

Description	Package Contains	Order Number
Regular	10	A00047

**Seating Lug/Hook**

- Serves as lingual hook or molar hook for cross elastics and as a band seating lug

Description	Package Contains	Order Number
Bicuspid Molar Base	100	A00767