


## Altitude SL® Brackets

RMO®'s Altitude SL® Brackets are the new solution to low-friction, passive self-ligation technology that delivers trusted results and treatment efficiency. Fast and efficient treatment is attained by effortless door activation and optimal sliding mechanics. Altitude SL® offers reliable strength in construction combined with a smooth rounded profile for enhanced patient comfort.



**MATERIAL:** 17-4PH Stainless Steel Pad  
**BASE:** 304 Stainless Steel 80-gauge Mesh  
**PRESCRIPTIONS:** Roth & MBT™  
**SLOT SIZE:** .018 & .022  
**TYPE:** Passive Self-Ligation

### FEATURES

Smooth, Rounded Profile  
 Ball End Hook  
 Optional Ligation  
 Torque-in-Base  
 Two-Piece Mesh Construction  
 Smooth Rounded Slot Lead-in  
 Patented Rail Design  
 Door Notch  
 Door Closed Gingivally  
 Base Palmer Notation  
 17-4PH Grade SS

### BENEFITS

Enhanced patient comfort  
 Added patient comfort through minimized cheek irritation  
 Ideal finishing and ligation  
 Optimal slot alignment reduces guesswork and chair time  
 Consistent and predictable bonding  
 Easy wire placement and insertion for quicker wire changes  
 Lower friction and binding on wire  
 Provides reliable, stable door retention  
 Quick and easy to open and close  
 No instrument required  
 Keeps door closed positively with chewing  
 Easy, permanent bracket identification  
 Molded for high strength and rigidity