

Ricketts® Classic Bioprogressive® Prescription - Mandibular Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018
Central/ Lateral	0°	0°	0°		R/L	M02620
Cuspid	3°	3°	0°		R	M02644
	3°	3°	0°		L	M02645
	3°	3°	0°	Center	R	M02646
	3°	3°	0°	Center	L	M02647
Optional High Profile	2°	2°	0°		R	M02624
	2°	2°	0°		L	M02625
	2°	2°	0°	Center	R	M02626
	2°	2°	0°	Center	L	M02627
1st Bicuspid	0°	0°	0°		R	M02628
Gingival Offset	0°	0°	0°		L	M02629
	0°	0°	0°	Center	R	M02658
	0°	0°	0°	Center	L	M02659
	0°	0°	0°		R	M02630
2nd Bicuspid	0°	0°	0°		L	M02631
	0°	0°	0°	Center	R	M02660
	0°	0°	0°	Center	L	M02661
	-14°	0°	0°		R	M02632
Gingival Offset	-14°	0°	0°		L	M02633
	-14°	0°	0°	Center	R	M02662
	-14°	0°	0°	Center	L	M02663
	-14°	0°	0°		R	M02634
Gingival Offset	-14°	0°	0°		L	M02635
	-14°	0°	0°	Center	R	M02664
	-14°	0°	0°	Center	L	M02665



Ricketts® Classic Bioprogressive® Prescription - Kits

Description	.018
Ricketts® Classic Synergy® without hook, 5-5	KM02610
Ricketts® Classic Synergy® hook on cuspid, 5-5	KM02611
Ricketts® Classic Synergy® hook on cuspid/bicuspid, 5-5	KM02612