

Roth Prescription - Maxillary Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central	12°	5°	0°		R	M05300	M05400
	12°	5°	0°		L	M05301	M05401
Lateral	8°	9°	0°		R	M05302	M05402
	8°	9°	0°		L	M05303	M05403
Cuspid	0°	11°	0°		R	M05304	M05404
	0°	11°	0°		L	M05305	M05405
	0°	11°	0°	Distal	R	M05306	M05406
	0°	11°	0°	Distal	L	M05307	M05407
1st/2nd Bicuspid	-7°	0°	0°		R/L	M05308	M05408
	-7°	0°	0°	Distal	R	M05310	M05410
	-7°	0°	0°	Distal	L	M05311	M05411
Gingival Offset	-7°	0°	0°		R/L	M05350	M05450
	-7°	0°	0°	Distal	R	M05352	M05452
	-7°	0°	0°	Distal	L	M05353	M05453



Roth Prescription - Mandibular Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central/	-1°	2°	0°		R	M05320	M05420
Lateral	-1°	2°	0°		L	M05321	M05421
Cuspid	-11°	5°	0°		R	M05324	M05424
	-11°	5°	0°		L	M05325	M05425
	-11°	5°	0°	Distal	R	M05326	M05426
	-11°	5°	0°	Distal	L	M05327	M05427
1st Bicuspid	-17°	0°	0°		R	M05328	M05428
	-17°	0°	0°		L	M05329	M05429
	-17°	0°	0°	Distal	R	M05330	M05430
	-17°	0°	0°	Distal	L	M05331	M05431
Gingival Offset	-17°	0°	0°		R	M05360	M05460
	-17°	0°	0°		L	M05361	M05461
	-17°	0°	0°	Distal	R	M05362	M05462
	-17°	0°	0°	Distal	L	M05363	M05463
2nd Bicuspid	-22°	0°	0°		R	M05332	M05432
	-22°	0°	0°		L	M05333	M05433
	-22°	0°	0°	Distal	R	M05334	M05434
	-22°	0°	0°	Distal	L	M05335	M05435
Gingival Offset	-22°	0°	0°		R	M05364	M05464
	-22°	0°	0°		L	M05365	M05465
	-22°	0°	0°	Distal	R	M05366	M05466
	-22°	0°	0°	Distal	L	M05367	M05467



Roth Prescription - Kits

Description	.018	.022
Roth FLI® Twin, without hook, 5-5	KM05300	KM05400
Roth FLI® Twin, hook on cuspid, 5-5	KM05301	KM05401
Roth FLI® Twin, hook on cuspid/bicuspid, 5-5	KM05302	KM05402
Roth FLI® Twin, hook on cuspid/gingival offset bicuspids, 5-5	KM05305	KM05405