


## Mini-Taurus® Brackets

Mini-Taurus® is the original MIM manufactured bracket. Setting standards since 1980, this bracket is smoothing out orthodontics for patients and doctors worldwide. Mini-Taurus features rounded contours and gentle hooks. This smoother, more comfortable design is then precisely manufactured with our metal injection molding (MIM) system.



**MATERIAL:** 316L Stainless Steel

**BASE:** 316L Stainless Steel 80-gauge Foil Mesh

**PRESCRIPTIONS:** Roth, Standard Edgewise, Alexander

**SLOT SIZE:** .018 & .022

### FEATURES

Rounded surfaces and edges

MIM, 80-Gauge Mesh Pad

Color quadrant & Permanent ID marking

Center scribe

### BENEFITS

Optimal patient comfort and aesthetically pleasing

Consistent, strong, predictable bonding

Quick and accurate identification

Offers precise bracket placement