

Roth Prescription - Maxillary Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	Base	.018	.022
Central	12°	5°	0°		R		M01000	M01100
	12°	5°	0°		L		M01001	M01101
Lateral	8°	9°	0°		R		M01002	M01102
	8°	9°	0°		L		M01003	M01103
	8°	9°	0°	Distal	R		M01016	M01116
	8°	9°	0°	Distal	L		M01017	M01117
Cuspid	0°	11°	4° MO		R		M01004	M01104
	0°	11°	4° MO		L		M01005	M01105
	0°	11°	4° MO	Distal	R		M01006	M01106
	0°	11°	4° MO	Distal	L		M01007	M01107
1st/2nd Bicuspid	-7°	0°	2° DO		R		M01008	M01108
	-7°	0°	2° DO		L		M01009	M01109
	-7°	0°	2° DO	Mesial	R		M01010	M01110
	-7°	0°	2° DO	Mesial	L		M01011	M01111
	-7°	0°	2° DO		R	On Band	A06408	A06508
	-7°	0°	2° DO		L	On Band	A06409	A06509
	-7°	0°	2° DO	Mesial	R	On Band	A06410	A06510
	-7°	0°	2° DO	Mesial	L	On Band	A06411	A06511
Gingival Offset	-7°	0°	2° DO		R		M01050	M01150
	-7°	0°	2° DO		L		M01051	M01151
	-7°	0°	2° DO	Mesial	R		M01052	M01152
	-7°	0°	2° DO	Mesial	L		M01053	M01153



Roth Prescription - Mandibular Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	Base	.018	.022
Central/	-1°	2°	0°		R		M01020	M01120
Lateral	-1°	2°	0°		L		M01021	M01121
Lateral	-1°	2°	0°	Distal	R		M01038	M01138
	-1°	2°	0°	Distal	L		M01039	M01139
Cuspid	-11°	5°	2° MO		R		M01024	M01124
	-11°	5°	2° MO		L		M01025	M01125
	-11°	5°	2° MO	Distal	R		M01026	M01126
	-11°	5°	2° MO	Distal	L		M01027	M01127
1st Bicuspid	-17°	0°	4° DO		R		M01028	M01128
	-17°	0°	4° DO		L		M01029	M01129
	-17°	0°	4° DO	Distal	R		M01030	M01130
	-17°	0°	4° DO	Distal	L		M01031	M01131
	-17°	0°	4° DO		R	On Band	A06428	A06528
	-17°	0°	4° DO		L	On Band	A06429	A06529
	-17°	0°	4° DO	Distal	R	On Band	A06430	A06530
Gingival Offset	-17°	0°	4° DO	Distal	L	On Band	A06431	A06531
	-17°	0°	4° DO		R		M01058	M01158
	-17°	0°	4° DO		L		M01059	M01159
	-17°	0°	4° DO	Distal	R		M01060	M01160
	-17°	0°	4° DO	Distal	L		M01061	M01161



Continued on next page