


## VersaT® Brackets

The patented VersaT® brings together all the features for reduced friction in a twin design along with a vertical slot for auxiliaries. VersaT®, by using the vertical slot, allows the clinician to adjust the tip, torque, and rotation and manage space to fit the individualized treatment plan for each patient.



**MATERIAL:** 316L Stainless Steel

**BASE:** 316L Stainless Steel 80-gauge Foil Mesh

**PRESCRIPTIONS:** Roth

**SLOT SIZE:** .018 & .022

### FEATURES

Rounded surfaces and edges

MIM, 80-Gauge Mesh Pad

Color quadrant & Permanent ID marking

Center scribe

### BENEFITS

Optimal patient comfort and aesthetically pleasing

Consistent, strong, predictable bonding

Quick and accurate identification

Offers precise bracket placement