

Roth Prescription - Maxillary Brackets

| Tooth | Torque | Angulation | Offset | Hook | R/L | .018 | .022 |
|----------|--------|------------|--------|--------|-----|--------|--------|
| Central | 12° | 5° | 0° | | R | Q03010 | Q03030 |
| | 12° | 5° | 0° | | L | Q03011 | Q03031 |
| Lateral | 8° | 9° | 0° | | R | Q03012 | Q03032 |
| | 8° | 9° | 0° | | L | Q03013 | Q03033 |
| Cuspid | 0° | 11° | 4° MO | | R | Q03018 | Q03038 |
| | 0° | 11° | 4° MO | | L | Q03017 | Q03037 |
| | 0° | 11° | 4° MO | Distal | R | Q03014 | Q03034 |
| | 0° | 11° | 4° MO | Distal | L | Q03015 | Q03035 |
| 1st/2nd | -7° | 0° | 0° | | R/L | Q03019 | Q03039 |
| Bicuspid | -7° | 0° | 0° | Center | R/L | Q03016 | Q03036 |



Roth Prescription - Mandibular Brackets

| Tooth | Torque | Angulation | Offset | Hook | R/L | .018 | .022 |
|---------------------|--------|------------|--------|--------|-----|--------|--------|
| Central/ Lateral | 0° | 0° | 0° | | R/L | Q03029 | Q03049 |
| Cuspid | -11° | 5° | 2° MO | | R | Q03025 | Q03045 |
| | -11° | 5° | 2° MO | | L | Q03026 | Q03046 |
| | -11° | 5° | 2° MO | Distal | R | Q03021 | Q03041 |
| | -11° | 5° | 2° MO | Distal | L | Q03022 | Q03042 |
| 1st | -17° | 0° | 0° | | R/L | Q03027 | Q03047 |
| Bicuspid | -17° | 0° | 0° | Center | R/L | Q03023 | Q03043 |
| 2nd | -22° | 0° | 0° | | R/L | Q03028 | Q03048 |
| Bicuspid | -22° | 0° | 0° | Center | R/L | Q03024 | Q03044 |



Roth Prescription - Kits

| Description | .018 | .022 |
|---|---------|---------|
| Roth Singature III hook on cuspid, 3-3 U | KQ03013 | KQ03033 |
| Roth Singature III without hook, 3-3 U/L | KQ03020 | KQ03040 |
| Roth Singature III hook on cuspid, 3-3 U/L | KQ03022 | KQ03042 |
| Roth Singature III without hook, 5-5 U/L | KQ03024 | KQ03044 |
| Roth Singature III hook on cuspid, 5-5 U/L | KQ03025 | KQ03045 |
| Roth Singature III hook on cuspid/bicuspid, 5-5 U/L | KQ03026 | KQ03046 |