

tads

ACCESSORIES: EYELETS, HOOKS, PLIERS

**TAD NiTi Coil Springs with Eyelets**

- Stainless steel eyelets are attached to each end of the closed NiTi coil spring to easily engage various appliances
- Wire diameter of .010" and inside coil diameter of .030"

Description	Package Contains	Order Number
7mm	10	G00304
9mm	10	G00301
12mm	10	G00302
15mm	10	G00303

**TAD NiTi Coil Springs with Oval Eyelets**

- Stainless steel eyelets attached to NiTi Coil Springs
- Pear-shaped opening allows for easy and secure attachment to Button Top TAD
- Wire diameter of .010" and inside coil diameter of .030"

Description	Package Contains	Order Number
8mm	10	G00305
10mm	10	G00306
12mm	10	G00307



A00130

A00129

Crimpable Hooks

- Designed for rectangular wire .016" x .022" through .019" x .025"
- Angled head improves access to the posterior region

Description	Package Contains	Order Number
Crimpable Hook	10	A00129
Long Crimpable Hook	10	A00130



Schweickhardt



Apex Series™

Crimpable Hook Pliers

- Helps to attach crimpable hooks securely to the arch wire
- Angled head improves access to the posterior region

Description	Package Contains	Order Number
Schweickhardt Crimpable Hook Plier	1 each	i00129
Apex Series™ Crimpable Hook Plier	1 each	T00130