




Instrument Rack


Package Contains	Order Number
1	I00042



Cephalometric Tracing Template

- Designed to provide the most anatomically correct presentation of typical teeth
- Includes: prediction table, frontal and lateral views of most teeth, pterygoid and condylar tracing shapes, metric scales, protractor and relationship guides
- Fits approximately 80% of the teeth without compensation, yet can also adjust to fit individual teeth.

Package Contains	Order Number
1 pkg	I00512



Penta-Morphic® Arch Selector

- Made of 8" x 10" Vinylite

Package Contains	Order Number
1	I00509



3D® Lingual Arch Selector

Package Contains	Order Number
1	I00505



3D® Bimetric Arch Selector

Package Contains	Order Number
1	I00506



3D® Maxillary Bimetric Distalizing Arch Selector

Package Contains	Order Number
1	I00507



Distalizing Arch Selector

Package Contains	Order Number
1	G00360



RMO® Instrument Lubricant

- Technical grade white oil that can be applied after drying and prior to sterilization
- Improves the instrument life
- Fast and economical solution that prevents rust and protects cutting edges
- Regular treatment reduces the chance of oxidation

Package Contains	Order Number
1	J00201