

instruments


APEX SERIES: WIRE BENDING PLIERS



Lingual Arch Forming Plier

- Plier makes consistent and repeatable bends
- Forms double-back and triple-back bends in .030" and .036" wire for use with lingual sheaths
- Jaws have four grooves and are sized for any type of larger wires
- Box joint

Wire Range	Package Contains	Order Number
.020" - .036" (0.508mm - 0.914mm) Round	1	T00420
.018" x .025" (0.457mm x 0.635mm) Rectangular		



Rectangular Arch Bending Plier (Tweed Style)

- Ideal for forming square and rectangular arch wires
- Thin jaws are parallel at .020" opening
- Forms offsets and 90° bends on all types of arch wire
- Box joint

Description	Package Contains	Order Number
.022" x .028" (0.559mm x 0.711mm) Rectangular	1	T00442



Crimpable Hook Plier

- Angled beak of the plier allows for precise placement of the hooks on the arch wire
- Hook is secure after on squeeze


Package Contains	Order Number
1	T00130



Step Plier

- Step plier is equipped with a double-sided offset beak for bayonet style step bend
- Creates right hand or left hand step-up or step-downs
- Box joint

Wire Range	Package Contains	Order Number
.012" - .020" (0.305mm - 0.508mm) Round	1	T00132
.022" x .028" (0.559mm x 0.711mm) Rectangular		
.012" - .020" (0.305mm - 0.508mm) Round	1	T00134
.022" x .028" (0.559mm x 0.711mm) Rectangular		



NiTi Crimping Plier

- Forms a smooth bend or crimp without nicking Nickel Titanium wire

Description	Package Contains	Order Number
.022" x .028" (0.559mm x 0.711mm) Rectangular	1	T00160