


### Disposable Impression Trays

- Made from FDA-approved styrene
- Features smooth, rounded edges, anterior retention slots and vertical perforations for enhanced material retention
- Anatomically correct design, with a palatal ridge for optimum fit, makes for greater patient comfort
- Ideal for use with RMO® Ortho-Jel® Alginate or any other impression material

Size	Color	Package Contains	Order Number: Maxillary
Small	Blue	10	J01233
Medium	Green	10	J01234
Large	Yellow	10	J01235
Assortment, Mixed		30 trays	K00455

Size	Color	Package Contains	Order Number: Mandibular
Small	Blue	10	J01243
Medium	Green	10	J01244
Large	Yellow	10	J01245
Assortment, Mixed		30 trays	K00455



### Aluminum Impression Trays

- Bends easily for a custom fit
- Smooth edges for optimal patient comfort
- Sterilizable

#### Maxillary

Size	Color	Package Contains	Order Number Maxillary
1 (57mm x 48mm)	Pink	1	J00304
2 (62mm x 54mm)	Blue	1	J00305
3 (68mm x 60mm)	Red	1	J00306
4 (75mm x 65mm)	Gold	1	J00307
5 (84mm x 70mm)	Violet	1	J00308

#### Mandibular

Size	Color	Package Contains	Order Number Mandibular
1 (62mm x 44mm)	Pink	1	J00364
2 (70mm x 48mm)	Blue	1	J00365
3 (75mm x 50mm)	Red	1	J00366
4 (82mm x 53mm)	Gold	1	J00367
5 (88mm x 55mm)	Violet	1	J00368



Description	Package Contains	Order Number
Aluminum impression tray kit	assortment of 10, one of each max. and man. sizes	K00566

# intraoral & extraoral


EXTRAORAL PRODUCTS: ALGINATE, CHEEK RETRACTORS, MOUTH GUARDS, RETAINER CASE



## Ortho-Jel® Alginate

- Fast Setting
- Vanilla mint scent
- 3 phase color indicator

Description	Package Contains	Order Number
Ortho-Jel®	(2) 450 gram bags	J00250
	(24) 450 gram bags	J00251



## Cheek Retractors

- Suction rings allow for use of hands-free suction for increased efficiency
- The addition of handles creates on-demand ability to increase clearance in hard to reach posterior areas
- Two sizes create additional sizing options for patient comfort
- Sterilization and autoclave-approved to increase office efficiency
- Wide opening creates a full view for bonding of anterior and posterior teeth

Description	Package Contains	Order Number
Cheek Retractor with Handle	1	J00405
Cheek Retractor, Adult	1	J00406
Cheek Retractor, Child	1	J00407



## Mouth Guards

- Comes in standard black and 3 fun colors/flavors
- Flavored Mouth Guards designed to fit with braces
- Safety design creates spacing between jaw and skull for protection from injuries
- Quick-release tether for flexibility of use across several sports
- Comes with reusable storage case for mobility
- Custom molding can be done using boil-and-bite procedure (standard black only)

Description	Color	Package Contains	Order Number
Standard	Black	1	J00555
Braces, Blue Raspberry	Blue	1	J00556
Braces, Fruit Punch	Red	1	J00557
Braces, Bubblegum	Pink	1	J00559



## Retainer Case

- Made of high impact polypropylene, prevents breakage
- Easy-open finger holds with durable hinge and positive-click closure
- Built-in keyhole to carry on key chain or backpack
- Ventilation holes that prevent odors
- Standard size holds upper or lower retainers or a positioner
- Lab size holds bionators or special lab appliances

Description	Color	Package Contains	Order Number
Standard Size	Tropical Pink	20	J01320
Standard Size	Melon	20	J01321
Standard Size	Metallic Black	20	J01323
Lab Size	Berry Blue	12	J01324