









1st Molar - Mandibular Single Tubes

	Description	Style	Hook	R/L	.018	.022
	Standard Edgewise, convertible	Weld		R/L	A04840	A04850
	Standard Edgewise, convertible	Bond		R/L	A08040	A08050
	Standard Edgewise, convertible	Weld	Mesial	R	A04842	A04852
	Standard Edgewise, convertible	Weld	Mesial	L	A04843	A04853
	Standard Edgewise, convertible	Bond	Mesial	R	A08042	A08052
	Standard Edgewise, convertible	Bond	Mesial	L	A08043	A08053
 FLI	FLI® Standard Edgewise	Weld	Mesial	R	A06778	A06782
	FLI® Standard Edgewise	Weld	Mesial	L	A06779	A06783
	FLI® Standard Edgewise	Bond	Mesial	R	A08778	A08782
	FLI® Standard Edgewise	Bond	Mesial	L	A08779	A08783
 FLI	FLI® Tube	Weld	Mesial	R	A06766	A06768
	FLI® Tube	Weld	Mesial	L	A06767	A06769
	FLI® Tube	Bond	Mesial	R	A08766	A08768
	FLI® Tube	Bond	Mesial	L	A08767	A08769
 FLI	FLI® MBT	Weld	Mesial	R	A05426	A05526
	FLI® MBT	Weld	Mesial	L	A05427	A05527
	FLI® MBT	Bond	Mesial	R	A08572	A08576
	FLI® MBT	Bond	Mesial	L	A08573	A08577
	SWLF®, convertible	Weld	Mesial	R	A05714	A05734
	SWLF®, convertible	Weld	Mesial	L	A05715	A05735
	SWLF®, convertible	Bond	Mesial	R	A08538	A08558
	SWLF®, convertible	Bond	Mesial	L	A08539	A08559
	Alexander, convertible not preangled	Weld	Mesial	R	A04384	A04394
	Alexander, convertible not preangled	Weld	Mesial	L	A04385	A04395
	Alexander, convertible not preangled	Bond	Mesial	R	A08184	A08194
	Alexander, convertible not preangled	Bond	Mesial	L	A08185	A08195
	Ricketts Bioprogressive®, without cap	Weld	Mesial	R	A04968	
	Ricketts Bioprogressive®, without cap	Weld	Mesial	L	A04969	
	Ricketts Bioprogressive®, without cap	Bond	Mesial	R	A08368	
	Ricketts Bioprogressive®, without cap	Bond	Mesial	L	A08369	

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-20°	0°	0°							
-20°	0°	0°							
-20°	0°	0°							
-20°	0°	0°							
-20°	0°	8° D							
-20°	0°	8° D							
-20°	0°	8° D							
-20°	0°	8° D							
-22°	-6°	5° D							
-22°	-6°	5° D							
-22°	-6°	5° D							
-22°	-6°	5° D							
-24°	-5°	12° D							
-24°	-5°	12° D							
-24°	-5°	12° D							
-24°	-5°	12° D							

1st Molar - Mandibular Single Tubes

Description	Style	Hook	R/L	.018	.022
 Ricketts Bioprogressive®, convertible	Weld	Mesial	R	A04978	
Ricketts Bioprogressive®, convertible	Weld	Mesial	L	A04979	
Ricketts Bioprogressive®, convertible	Bond	Mesial	R	A08378	
Ricketts Bioprogressive®, convertible	Bond	Mesial	L	A08379	
 FLI® Tube	Weld	Mesial	R	A06740	A06744
FLI® Tube	Weld	Mesial	L	A06741	A06745
FLI® Tube	Bond	Mesial	R	A08740	A08744
FLI® Tube	Bond	Mesial	L	A08741	A08745
 FLI® One Piece	Bond	Mesial	R	A18758	A18762
FLI® One Piece	Bond	Mesial	L	A18759	A18763
 Bench Bioprogressive®, without cap	Weld	Mesial	R	A01196	
Bench Bioprogressive®, without cap	Weld	Mesial	L	A01197	
Bench Bioprogressive®, without cap	Bond	Mesial	R	A08296	
Bench Bioprogressive®, without cap	Bond	Mesial	L	A08297	
 Bench Bioprogressive®, convertible	Weld	Mesial	R	A01198	
Bench Bioprogressive®, convertible	Weld	Mesial	L	A01199	
Bench Bioprogressive®, convertible	Bond	Mesial	R	A08298	
Bench Bioprogressive®, convertible	Bond	Mesial	L	A08299	
 Roth, convertible, not preangled	Weld	Mesial	R	A05666	A05694
Roth, convertible, not preangled	Weld	Mesial	L	A05667	A05695
Roth, convertible, not preangled	Bond	Mesial	R	A08638	A08640
Roth, convertible, not preangled	Bond	Mesial	L	A08639	A08641
 Andrews, convertible, not preangled	Weld	Mesial	R	A06120	A06160
Andrews, convertible, not preangled	Weld	Mesial	L	A06121	A06161
Andrews, convertible, not preangled	Bond	Mesial	R	A08420	A08460
Andrews, convertible, not preangled	Bond	Mesial	L	A08421	A08461
 FLI Roth, not preangled	Weld	Mesial	R	A05608	A05610
FLI Roth, not preangled	Weld	Mesial	L	A05609	A05611
FLI Roth, not preangled	Bond	Mesial	R	A08208	A08210
FLI Roth, not preangled	Bond	Mesial	L	A08209	A08211

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-24°	-5°	12° D							
-24°	-5°	12° D							
-24°	-5°	12° D							
-24°	-5°	12° D							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	5° D							
-25°	0°	5° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-30°	-1°	0°							
-30°	-1°	0°							
-30°	-1°	0°							
-30°	-1°	0°							
-30°	+2°	0°							
-30°	+2°	0°							
-30°	+2°	0°							
-30°	+2°	0°							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							