




## SELF-LIGATING MOLAR TUBES

**Self-ligating passive clip for retaining arch wire produces less friction and lower force mechanics throughout treatment**

- No need to retie converted 6th tube when 2nd molars are in use
- Easier wire insertion with convertible clip and funneled mesial slot entrance
- Robust clip with an efficient key hole for dependable performance
- Full slot-width clip coverage delivers exceptional rotational control
- Eliminates ligatures for increased lower sliding mechanics over traditional molar tubes
- Anatomical contoured base for better retention and unsurpassed performance
- Low-profile design provides maximum patient comfort
- Small and strong hooks for coil springs, chain elastics and other applications
- Time-saving benefits and ease of self-ligation
- Ultra-smooth steel construction with chromium cobalt clip
- Ease of use with color ID markings

**Easier wire insertion** - A funneled tube entrance reduces wire binding and allows for easier wire insertion

**Exceptional rotational control** - A full slot-width clip delivers exceptional rotational control




**Low-profile hooks** - Designed for patient comfort and improved hygiene

**Simplified clip opening** - Key hole simplifies opening the clip

# SELF-LIGATING MOLAR TUBES

## Another innovation from RMO!



A08840	SL MOLAR SGL MX6 R18 -10T 10DO HK	Black
A08841	SL MOLAR SGL MX6 L18 -10T 10DO HK	Green
A08842	SL MOLAR SGL MX6 R18 -14T 10DO HK	Black
A08843	SL MOLAR SGL MX6 L18 -14T 10DO HK	Green
A08844	SL MOLAR SGL MN6 R18 -10T HK	Red
A08845	SL MOLAR SGL MN6 L18 -10T HK	Yellow
A08846	SL MOLAR SGL MN6 R18 -20T HK	Red
A08847	SL MOLAR SGL MN6 L18 -20T HK	Yellow
A08848	SL MOLAR SGL MN6 R18 -25T HK	Red
A08849	SL MOLAR SGL MN6 L18 -25T HK	Yellow
A08850	SL MOLAR SGL MX6 R22 -10T 10DO HK	Black
A08851	SL MOLAR SGL MX6 L22 -10T 10DO HK	Green
A08852	SL MOLAR SGL MX6 R22 -14T 10DO HK	Black
A08853	SL MOLAR SGL MX6 L22 -14T 10DO HK	Green
A08854	SL MOLAR SGL MN6 R22 -10T HK	Red
A08855	SL MOLAR SGL MN6 L22 -10T HK	Yellow
A08856	SL MOLAR SGL MN6 R22 -20T HK	Red
A08857	SL MOLAR SGL MN6 L22 -20T HK	Yellow
A08858	SL MOLAR SGL MN6 R22 -25T HK	Red
A08859	SL MOLAR SGL MN6 L22 -25T HK	Yellow



To order contact customer service or call:  
**800.525.6375**