

bracket systems

MISCELLANEOUS BRACKETS: PERAHIA SANCHEZ, I-TTR



Perahia Sanchez Prescription - Maxillary/Mandibular Brackets

Tooth	Torque	Angulation	Vertical Slot	Offset	Style	R/L	.018
Central/ Lateral	0°	0°	.019 x .020	0°		R/L	M01777
Cuspid/ Bicuspid	0°	0°	.019 x .020	0°		R/L	M01779



i-TTR Lingual System

The intelligent lingual treatment with torque, tip and rotation (i-TTR) control. The bracket presents two slots for round arches and a rectangular central slot for an auxiliary insert, which permit a progressive modular control of all tooth movement. This can be used for all teeth, including molars.

Description	Package Contains	Order Number
i-TTR Lingual Bracket	10	M00700
i-TTR Lingual Bracket Kit (Double brackets for first and second molars)	36	KM00700



Light Wire Insert Pin

For A01777 Bracket. Specially-hardened stainless steel light wire insert pin designed for .016 (0.406mm) and .018 (0.457mm) wires. Pin provides light-wire functions during treatment. Also keeps cement from entering bandside slot.

A00739 PKG 100 Approx.