

adhesives

LIGHT AND PRIMER ACTIVATED ADHESIVES: SYRINGES, KITS



TruLock® Light Activated Adhesive

TruLock® Light Activated Adhesive is a single paste resin-based, fluoride containing bracket adhesive that requires no mixing, resulting in reduced waste, improved performance consistency and no working time constraints.

- Application: Bonds metal, ceramic and plastic/composite materials to enamel
- Shelf Life: 30 months from manufacture date
- Curing: Minimum output 600 mW/cm² and wavelength 400-500 nanometers

Description	Package Contains	Order Number
Light Activated Adhesive (5g Syringe)	4 each	J04000
Light Activated Adhesive (40 X .35g Capsule) 14g	40 each	J04001
Light Activated Bonding Resin 7ml	1 each	J04002
Description	Package Contains	Order Number
Light Activated Kit - 4 x 5g Syringes, 2 x 3ml Bonding Resin, 7ml Etchant Liquid, Brushes, Spatula, Mixing Pad	1 each	K00070
Light Activated Kit - 40 x 0.35g Capsules, 7ml Bonding Resin, 2 x 4g Etchant Gel Syringes, Tips, Brushes, Mixing Pad	1 each	K00071
Light Activated Mini Kit - 5g Syringe, 3ml Bonding Resin, 3ml Etchant Liquid, Brushes, Spatula, Mixing Pad	1 each	K00072
Light Activated Mini Kit - 5g Syringe, 3ml Bonding Resin, Brushes, Spatula, Mixing Pad, (NO ETCHANT)	1 each	K00073
Light Activated Kit - 4 x 5g Syringes, 2 x 3ml Bonding Resin, Brushes, Spatula, Mixing Pad, (NO ETCHANT)	1 each	K00074



TruLock® Primer Activated Adhesive

TruLock® Primer Activated Adhesive is a paste-primer orthodontic formula. The unique filler and primer combination produces an ideal bond and a smooth viscosity minimizing bracket drift or movement.

- Application: Bonds metal, ceramic, and plastic/composite materials to enamel
- Shelf Life: 30 months from manufacture date
- Curing: Self cured with paste and primer

Description	Package Contains	Order Number
Primer Activated Adhesive (jar), 14g	1 each	J04010
Primer Activated Adhesive (3.5g Syringe)	40 each	J04011
Primer Activated Adhesive (40 X .35g Capsule), 14g	1 each	J04012
Primer Activated Primer, 12ml	1 each	J04013
Description	Package Contains	Order Number
Primer Activated Kit - 14g Adhesive Jar, 12ml Primer, 7ml Etchant Liquid, Brushes, Spatulas, Mixing Pad	1 each	K00080
Primer Activated Kit - 4 x 3.5g Syringes, 12ml Primer, 7ml Etchant Liquid, Brushes, Spatulas, Mixing Pad	1 each	K00081
Primer Activated Kit - 40 x 0.35g Capsules, 12ml Primer, 7ml Etchant Liquid, Brushes, Mixing Pad	1 each	K00082
Primer Activated Mini Kit - 3.5g Syringe, 5ml Primer, 3ml Etchant Liquid, Brushes, Spatulas, Mixing Pad	1 each	K00083
Primer Activated Mini Kit - 3.5g Syringes, 5ml Primer, Brushes, Spatulas, Mixing Pad (NO ETCHANT)	1 each	K00084
Primer Activated Kit - 4 x 3.5g Syringes, 12ml Primer, Brushes, Spatulas, Mixing Pad, (NO ETCHANT)	1 each	K00085